



Specification Sheet

NovAseptic® Valve Inductive Position Indicator Unit

Catalogue Number PA1251/SP Single Position

Inductive Position Indicator Unit





Nominal Dimensions in mm (in.)


PA##	Α	В	С	D	PA ##/100	PA ##/110	PA ##/111
PA1251/SP	53.0 (2.087)	5.0 (0.197)	30.0 (1.181)	30.9 (1.217)	PA12/100	PA12/110	PA12/111
PA1251/SP	53.0 (2.087)	5.0 (0.197)	30.0 (1.181)	30.9 (1.217)	PA18/100	PA18/110	PA18/111
PA1251/SP	53.0 (2.087)	5.0 (0.197)	30.0 (1.181)	30.9 (1.217)	PA25/100	PA25/110	PA25/111
PA1251/SP	53.0 (2.087)	5.0 (0.197)	30.0 (1.181)	30.9 (1.217)	PA38/100	PA38/110	PA38/111
PA1251/SP	53.0 (2.087)	5.0 (0.197)	30.0 (1.181)	30.9 (1.217)	PA51/100	PA51/110	PA51/111

Net Weight (approximate)

Indicator	PA1251/SP
Weight kg (lb)	0.06 (0.13)

Material Body		Indicator	
	plastic in compliance with	plastic in compliance with	
Material Code	РОМ С		


Socket Connection: M12x1



Power Supply	24V DC		
Power Limits	10-30V DC		
Output Current	Max 100Ma in open or closed position, overload and short circuit protection		
No Load Current Consumption	<15mA		
Degree of Protection	IP65		

Example

Signal output as a function of the position of a normal open or normal closed pneumatic actuator when mounted according to the mounting instructions.



Specifications

Surface Roughness	Smooth	
Environment Temperature	0 to 40 °C (32 to 104 °F) Under limited time 70 °C (158 °F)	
Labeling	Each Indicator is labeled.	
Packaging	The Indicator is packaged in a plastic bag.	
Quality Control	Millipore Quality Assurance System guarantees the control and traceability of the product.	

FOR ADDITIONAL INFORMATION CALL YOUR NEAREST MILLIPORE OFFICE:

In the U.S. and Canada, call toll-free 1-800-MILLIPORE (1-800-645-5476)

In the U.S., Canada and Puerto Rico, fax orders to

1-800-MILLIFX (1-800-645-5439)



Millipore, NovAseptic and Advancing Life Science Together are registered trademarks of Millipore Corporation. The M mark is a trademark of Millipore Corporation.

Lit. No. SP1202EN00 Rev. A 8/10

©2010 Millipore Corporation, Billerica MA, 01821 USA. All rights reserved.